

J.T.Baker[®] brand



High-purity, cGMP-manufactured Isopropyl-beta-D-Thiogalactopyranoside (IPTG)



HIGHLIGHTS

- cGMP manufactured under ICHQ7
- High-purity (>98% by HPLC)
- Low microbial, trace metals and endotoxin content
- Consistent lot-to-lot protein induction performance
- TSE/BSE-free
- 1,4-Dioxane-free

Reduce risk in your
bioprocessing applications with
Avantor J.T.Baker brand cGMP
manufactured IPTG

J.T.Baker brand cGMP manufactured IPTG is performance tested to assure consistent lot-to-lot biological induction and processed to remove and control product/process related impurities, with minimal lot to lot variability for improved protein expression.

Packaged in moisture permeation resistant controlled packaging, J.T.Baker brand cGMP manufactured IPTG is fully validated and documented through Batch Process Records (BPR), including analytical test data, stability protocol and impurity profile.

cGMP manufactured IPTG facilitates constant and consistent recombinant protein production resulting in more predictable yields.

APPLICATION INFORMATION

IPTG is a biotech reagent for protein production, used during microbial fermentation to induce the expression of recombinant proteins controlled by lac operon.

As a molecular mimic of allolactose, IPTG is a non-hydrolyzable galactose derivative that induces transcription of the gene in lac operon by binding to and inactivating the lac repressor. Unlike lactose, IPTG is not metabolized by cells, thus eliminating the need for continuous monitoring of the inducer concentration during protein expression. Additionally, because IPTG concentration remains constant, high titer production and consistent protein expression of genes controlled by the lac operon can be achieved.

The optimum concentration of IPTG used for inducing protein expression should be empirically determined as it can depend upon multiple factors, such as the strength of induction required, cellular genotype, etc.

IPTG is also a commonly used inducer in conjunction with X-Gal in blue-white color screen for rapid identifying recombinant bacterial colonies.

Additionally, unlike lactose, IPTG is non-animal derived eliminating the need for TSE documentation/testing.

PRODUCT NUMBERS

Pk	Cat. No
10 g	V264-01
100 g	V264-02
1 kg	V264-05

For more information, contact your Avantor account manager today.



VWR.COM

AUSTRIA

VWR International GmbH
Tel.: +43 1 97 002 0
info.at@vwr.com

BELGIUM

VWR International bvba
Tel.: +32 (0) 16 385 011
vwr.be@vwr.com

CZECH REPUBLIC

VWR International s. r. o.
Tel.: +420 321 570 321
info.cz@vwr.com

DENMARK

VWR International A/S
Tel.: +45 43 86 87 88
info.dk@vwr.com

FINLAND

VWR International Oy
Tel.: +358 (0) 9 80 45 51
info.fi@vwr.com

FRANCE

VWR International S.A.S.
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00 (international)
info.fr@vwr.com

* 0,18 € TTC/min + prix appel

GERMANY

VWR International GmbH
Tel.: 0800 702 00 07* (national)
Tel.: +49 (0) 6151 3972 0 (international)
info.de@vwr.com

*Freecall

HUNGARY

VWR International Kft.
Tel.: +36 52 521130
info.hu@vwr.com

IRELAND / NORTHERN IRELAND

VWR International Ltd / VWR
International (Northern Ireland) Ltd
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY

VWR International S.r.l.
Tel.: +39 02 3320311
info.it@vwr.com

THE NETHERLANDS

VWR International B.V.
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

NORWAY

VWR International AS
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND

VWR International Sp. z o.o.
Tel.: +48 58 32 38 200
info.pl@vwr.com

PORTUGAL

VWR International - Material de
Laboratório, Lda
Tel.: +351 21 3600 770
info.pt@vwr.com

SPAIN

VWR International Eurolab S.L.
Tel.: +34 902 222 897
info.es@vwr.com

SWEDEN

VWR International AB
Tel.: +46 (0) 8 621 34 00
kundservice.se@vwr.com

SWITZERLAND

VWR International GmbH
Tel.: +41 (0) 44 745 13 13
info.ch@vwr.com

UK

VWR International Ltd
Tel.: +44 (0) 800 22 33 44
uksales@vwr.com

CHINA

VWR (Shanghai) Co., Ltd
Tel.: 400 821 8006
Email: info_china@vwr.com

INDIA

VWR Lab Products Private Limited
Tel.: +91 80 28078400
vwr_india@vwr.com

MIDDLE EAST & AFRIKA

VWR International FZ-LLC
Tel.: +971 4 5573271
info.me@vwr.com

SINGAPORE

VWR Singapore Pte Ltd
Tel.: +65 6505 0760
sales.sg@vwr.com